



REHASENSE



# BACKUP

## SHOWERCHAIRS

Backup – a powder coated, height adjustable aluminium showerchair with fixed handles. Available either with a fixed or a knock-down frame with protruding armrests for better grip.

# BACKUP

Backup is CE-marked and fulfills the requirements of EU directive 93/42 EEC for medical devices.



REHASense

## SHOWERCHAIRS

### HEIGHT ADJUSTABLE

Seat height 44–64 cm, can be adjusted to 8 different positions with 2,5 cm steps by the use of snaplock clips, so it is suitable also for tall users. Seat width and depth: 40 x 40 cm.

### REMOVABLE BACKREST

Removable and length adjustable backrest is an accessory on the fixed frame showerchairs and not installed as standard. With the knock-down version the back comes as standard with all the parts in a one box. The backrest can be adjusted to 2 different positions in increments of 1 cm. Chairs can be used with or without a backrest.

### DURABLE & STRONG

The construction has been tested and approved on users with a maximum weight of 150 kg.



Fixed frame: 3 models, without backrest as standard.



### KNOCK-DOWN FRAME

Chairs are available either as a fixed version or as a knock-down version. The knock-down version comes with protruding armrests for better grip and always with the backrest.

### TECHNICAL INFORMATION

Max. user weight: 150 kg

Colours:

Frames: White or grey

Seat & back: Grey or black

Recommended user height: 130–210 cm

Width x length x height: without backrest:

56 x 53 x 60–80 cm

with backrest:

56 x 59–60 x 80–100 cm

Width between armrests: 42 cm

Seat height:

44–64 cm (every 2,5 cm)

Total weight: 4,8 kg

(incl. backrest at 2,6 kg)

Material:

Non-corrosive, powder coated aluminium & PU (polyurethane)



Knock-down version has small delivery dimensions 50 x 45 x 40 cm.



With one fine tuning leg for final, precise adjustment and levelling of the chair.



All legs with 8 x 2,5 cm increments to adjust seat height between 44–64 cm. Rubber ferrules prevents the chair from slipping on a slippery surface

Configurations item no.:

Frame	Armrest type	Stackable	Seat			# in box
			Standard		Hygienic	
			grey	black	grey	
Fixed	non-protruding	stackable	BUWB3001	BUWB3002	BUWB3001H	4
Knock-down	protruding	non-stackable	BUWB3111	BUWB3112	BUWB3111H	1

Back			# in box
grey	black		
BFSC3001	BFSC3002	optional	4
		included	-



REHASense  
Authorised Reseller

Copyright Reha-Care-Pro © 2018